

Supplementary Information

Development and application of a high-throughput method for the purification and analysis of surface carbohydrates from *Klebsiella pneumoniae*

Francesca Nonne¹, Mariagrazia Molfetta¹, Rebecca Nappini^{1,2}, Chiara La Guidara¹, Roberta Di Benedetto¹, Siwaphiwe Mfana³, Barbara Bellich^{4*}, Maria Michelina Raso¹, Gianmarco Gasperini⁵, Renzo Alfini¹, Paola Cescutti², Francesco Berlanda Scorza¹, Neil Ravenscroft³, Francesca Micoli¹, Carlo Giannelli¹

¹GSK Vaccines Institute for Global Health, Siena, Italy

²Department of Life Science, University of Trieste, Trieste, Italy

³Department of Chemistry, University of Cape Town, South Africa.

⁴Department of Advanced Translational Diagnostics, Institute for Maternal and Child Health, IRCCS "Burlo Garofolo", Trieste, Italy

⁵GSK, Siena, Italy

* When the experiments were performed BB was working at the Department of Life Sciences, University of Trieste, Trieste, Italy

Table S1. Evaluation of OAc level after KAg purification

KAg type	Strain	Extraction method	OAc %
K2	NCTC13810	H ₂ O	~55%
	NCTC13810	AcOH	~45%
	NCTC11228	AcOH	~98%
	5765B	H ₂ O	~77%

Table S2: Complete ¹H and ¹³C NMR chemical shifts data of the sugar residues in the polysaccharide produced by strain 4998 (identified as O-acetylated K2) with glycosylation shifts shown in brackets below the carbon chemical shifts; the linkage carbons are underlined.

Residue	H1 C1	H2 C2	H3 C3	H4 C4	H5 C5	H6 C6	OAc
A →4)-α-D-Glcp-(1→	5.28 99.79 (6.80)	3.61 72.27 (-0.20)	3.86 72.39 (-1.39)	3.66 <u>80.13</u> (9.42)	4.05 71.21 (-1.16)	3.75; 3.79 60.87 (-0.97)	-
B α-D-GlcpA-(1→	5.17 102.04 (9.08)	3.53 72.65 (0.39)	3.79 74.40 (0.88)	3.51 72.72 (-0.19)	4.10 73.45 (0.98)	- 176.89 (-0.53)	-
C →3,4)-β-D-Manp6Ac-(1→	4.78 100.86 (6.31)	4.23 71.42 (-0.71)	3.86 <u>81.14</u> (7.11)	4.12 <u>73.21</u> (5.52)	3.80 73.79 (-3.21)	4.66; 4.31 63.51 (1.52)	2.17 21.13 174.61
D →3)-β-D-Glcp-(1→	4.47 103.20 (6.36)	3.39 73.03 (-2.17)	3.63 <u>83.92</u> (7.16)	3.56 70.91 (0.20)	3.50 77.18 (0.42)	3.75; 3.97 61.76 (-0.08)	-




Figure S1. Western Blot analyses shows recognition of K62 and K64 samples by an anti-MrkA antibody confirming the presence of MrkA protein in the purified KAgS.




Figure S2. HPLC-SEC profiles highlight differences in the size of purified K-Ag obtained from bacterial extraction with and without AcOH.




Figure S3. Impact of AcOH on K2 (A) and K64 (B) PS structure integrity by NMR analysis





Figure S4. ^1H - ^1H NMR COSY (red)/TOCSY (black) overlay of the anomeric region of K-Ag from strain 4998 recorded at 600 MHz and 343 K. The crosspeaks for the O-acetylated tetrasaccharide RU are labeled (A = α -Glc, B = α -GlcA, C = β -Man6Ac and D = β -Glc).






Figure S7. The HSQC-DEPT spectrum of O-Ag from strain 7008B recorded at 600 MHz and 323 K, methylene signals are inverted (in red). All the pentasaccharide repeat unit proton/carbon crosspeaks have been labeled according to the carbon atom of the corresponding residue in the RU (A = β -3Gal f , B = α -3Gal p , C = α -3,4Gal p , D = α -Gal p and E = β -3Gal p). Additional peaks are due to buffer and glycerol (Gro).